

We recommend this product is installed by a fully qualified tradesperson.

Tools Required (not supplied)

Pencil; Tape measure; Screwdriver appropriate to fixings; Spirit level; Drill; Masonry drill bit appropriate to fixings, Ceramic drill bit (optional) appropriate to fixings.

Home Safety Advice

• It is recommended that two people fit this product to ensure that it is fitted safely.

• Always take care when using power tools, particularly in the bathroom.

• Always check for hidden cables and pipework before drilling and take extreme care if there is any water in the working area.

• It is advisable to use a residual current circuit breaker (RCCB).

• Always wear suitable eye protection when drilling.

• If the product is to be fitted on a ceramic tiled wall, a ceramic drill bit should be used.

• Always ensure that the drill hole passes through the central tiled area rather than through the grouted area. To prevent damage to the tile, mask the area around the hole with tape before drilling.

• Do not strike the product with hard or sharp objects.

• Warning! This product has been designed as a seat only and should not be used for any other purposes.

Maximum Weight - 160 kg evenly distributed about the seat.

Product Location

Shower seats are only as strong as the wall to which they are fixed and the fixings that are used. For different wall types you will need to ensure that you are using the correct type of fixing for the material of the wall. The addition of backboards on the outside of the wall or appropriate strengthening of internal structures may be required to provide adequate wall strength. Product Care

To retain the best quality finish, clean product regularly with a soft damp cloth. Do not use abrasive or chemical cleaners as these will damage the product. Dry the shower seat after use.

Note: If the wall is made of lightweight thermal blocks, aerated hollow brickwork, or prefabricated partitions, even the most appropriate fixings may not be able to withstand the loads which can be suddenly applied to shower seats.

Recycling and disposal

Please recycle the packaging in accordance with your local government regulations on waste handling and follow the same advice at the end of the products' life.

For more Information

Visit: victorianplumbing.co.uk

PARTS LIST		
ITEM	QTY	PART NUMBER
1	1	SEAT ASSEMBLY
2	8	6x60 WALL SCREW
3	8	GREY WALL PLUG

INSTRUCTION SUITABLE FOR SSW001



The fixing wall that you intend to mount the seat onto must be a solid wall. Alternatively, you can use a partition wall with suitable fixing points, such as a noggin or baton.

The fixing wall must have a fully sealed and cured covering such as tiles.

Take appropriate health and safety precautions when installing the product.







victorian plumbing.co.uk

3

1

To ensure the correct installation, it is recommended using a qualified, competent installer to install the product.

The shower tray and wall coverings, such as tiles, should be fully installed, sealed and cured.







2

Place the seat at the desired height and in the desired location. The legs are adjustable by unclipping the pin, adjusting the leg and then re-inserting the pin, In addition, one of the feet on the lower legs is height adjustable via a screw thread.



Make sure that the seat is level and stable in all planes in its lowered position. Mark the holes onto the wall.





Custom<1 Care Guide

3

- It is important to clean your shower seat after every use. wipe down all plastic and metal Parts with a soft cloth.
- Do not allow lime scale to build up in hard water areas.
- Ensure that any residual water is removed after every use.
- Do not use abrasive products on any surfaces.
- Bathroom products are vulnerable to add attack and some household consumables such hair dye, disinfectants, denture cleaners etc can contain such adds. Contact with these materials should be avoided.
- This instruction guide is for illustration purposes and may not be an exact representation of the physical product.
- For full terms and conditions, please refer to the manufacturer